

»» My new NUMBER ONE  
for monophasic impressions ««

# VSXE® ONE




Changed  
Cartridge System



SOOO RELIABLY TRANSFERRED

PRECISION IMPRESSIONS



KETTENBACHDENTAL  
Simply intelligent

# VSXE® ONE – UNIQUE MATERIAL COMBINATION

## VSXE® One:


Start now with VSXE® One, specifically designed for **monophase** impressions, such as **implant** impressions.

What sets VSXE® One apart, aside from its unique name, is its **explicit focus on fixation impressions**, with its product properties precisely tailored to this application.

## Indications / Techniques

VSXE® One is optimized for:

- Implant impressions
- Fixation impressions
- Functional impressions



## VSXE® ONE, THE NEW ULTIMATE FORMULA:



### VIGOROUS

Secure transfer due to **exceptionally high final hardness (shore-A 65)**, ensuring the reliable fixation of implant posts.



### SEAMLESS

Excellent flowability resulting from the formulation's unique combination providing **optimal coverage** of transfer posts.



### X-FORMULA

**Extended working time** (2:30 min.) and **quick setting time** (5:15 min.), preventing potential deformations during the setting phase and promoting **patient comfort**.



### EFFORTLESS

**Easy removal** from the model and the patient's mouth, safeguarding against unintended extractions



ADVANTAGES  
THAT **STICK**

## = Your new #ONE for monophase impressions!

- ✓ **Strong thixotropy** for an ideal pressure build-up
- ✓ **Stable in storage over a long period of time** for repeated pouring of the impression
- ✓ **Tasteless and odorless** for the well-being of the patients
- ✓ **No workflow change** in the practice or in the laboratory!

# VSXE® ONE – THE RANGE



## VSXE® One

380 ml **Intro pack** 5:1  
380 ml cartridge, 10 ml adhesive,  
10 dynamic mixers, 1 application syringe

REF 14732

760 ml **Refill pack** 5:1  
2 × 380 ml cartridges

REF 14733

## Accessories

**Universal Adhesive**  
10 ml bottle

REF 14206

**Dynamic Mixers**  
45 pieces

REF 17900

**Sympress Dispenser**  
Electrical dosing and  
mixing device

REF 35910

**Applyfix 5**  
2 application syringes,  
12 syringe tips

REF 17204



Made  
in  
Germany

## NOTES

---



---



---



---



---

Kettenbach SNC  
7, Place de la Gare  
57200 Sarreguemines · France  
info@kettenbach.com  
www.kettenbach-dental.com



1410734\_3524

**KETTENBACHDENTAL**  
Simply intelligent